

Experiment 6: Synthesis of 3-(4-Methoxyphenyl)pyridine using $\text{Pd}(\text{Ph}_3)_4$

Jessica Persaud, Jenna Pralat, Shania Taylor, Max Shi
We pledge our honor we have abided by the Stevens Honor Code.

Procedure


Preparation of Pyridylzincbromide reagent

1. Charge a 250 mL RBF with active zinc (9.81 g) in 100 mL THF
2. Add 3-Bromopyridine (15.8 g) to RBF via cannula while being stirred at room temperature
3. Stir mixture at refluxing temperature and monitor the oxidative addition by GC analysis
4. Let settle overnight and then transfer supernatant to 500 mL bottle
5. Dilute with 200 mL of THF

Coupling of Pyridylzincbromide with 4-iodoanisole

1. Add 0.10 g of $\text{Pd}[\text{P}(\text{Ph}_3)]_4$ to a 50 mL RBF
2. Add 20 mL of solution from part A and 1.80 g of 4-Iodoanisole dissolved in 10 mL THF
3. Stir at room temp for 1 hour then quench with saturated NH_4Cl solution
4. Extract with 30 mL of ethyl acetate, 3 times
5. Wash organic layers with saturated NaHCO_3 and brine
6. Dry with anhydrous MgSO_4
7. Use flash column chromatography (20%EtOAc/80% heptane)

Flow Chart




Preparation of Pyridyl Zinc Bromide reagent

1. Charge a 250 mL round bottom flask with 9.81 g zinc in 100 mL THF
2. Add 15.8 g 3-bromopyridine via cannula while being stirred at room temp.




Zinc




100 mL

9.81 grams



Stir mixture

15.8 grams




3-bromopyridine


Preparation of Pyridyl Zinc Bromide reagent

3. After addition, continue to stir at refluxing temperature, and monitor by GC analysis.

4. After being settled down overnight, transfer supernatant to 500 mL bottle and dilute with fresh THF to 200 mL.




Coupling of Pyridylzincbromide with 4-iodoanisole




Coupling of Pyridylzincbromide with 4-iodoanisole

Extract with 30 mL of ethyl acetate, 3 times



Organic layer is top

Wash with NaHCO_3 and brine




Add MgSO_4



Organic layer

Coupling of Pyridylzincbromide with 4-iodoanisole



Conclusion

Accomplished: In this experiment 3-(4-methoxyphenyl)pyridine was synthesized from 3-bromopyridine, zinc, and 4-iodoanisole using a Negishi coupling reaction with the $\text{Pd}(\text{Ph}_3)_4$ catalyst.

Learned: Through this experiment we learned how to synthesize 3-(4-methoxyphenyl)pyridine using the Negishi coupling reaction. We also learned how to use cannulation to transfer reagents between flasks under an inert atmosphere.

Issues: Potential issues can arise if cannulation was not performed properly. $\text{Pd}(\text{Ph}_3)_4$ is highly sensitive to air, thus any errors within its usage can result in a skewed product.

Recommendations / Applications: The $\text{Pd}(\text{Ph}_3)_4$ catalyst is commonly used in Suzuki, Negishi, and Heck coupling reactions. The Suzuki coupling reaction in particular is a scalable and cost effective process to produce intermediates for pharmaceuticals and fine chemicals. The coupling reaction requires common boronic acids which are less toxic and safer to the environment.

Post Lab Questions

- 1) List one advantage of using Arylzinc compounds compared to aryllithium or Grignard reagents for coupling reactions.
 - a) One advantage of using aryl zinc compounds over aryl lithium or grignard reagents is that they are less reactive, this is useful when partial reactions are desired, such as partial alkylation.
- 2) Give an example of another palladium containing catalyst and reaction involving that catalyst.
 - a) Another palladium catalyst is palladium(II) acetate which can be used for the Heck Reaction, as

