

# CH-412-LA Experiment 3

Oxidative and reductive properties of  $\text{H}_2\text{O}_2$  and  $\text{H}_2\text{S}$ ; precipitation and dissolution of sulfide compounds.

Max Shi, Brynn Froehlich, Erika McCarthy, Emily McGuire

*I pledge my honor that I have abided by the Stevens Honor System.*

# Procedure

## Oxidative and Reductive Activity of H<sub>2</sub>O<sub>2</sub>

1. Mix equal parts 0.2M Pb(NO<sub>3</sub>)<sub>2</sub> and H<sub>2</sub>S solution. After removing the supernatant, add 30% H<sub>2</sub>O<sub>2</sub> solution.
2. Add several drops of 0.01M KMnO<sub>4</sub> to a test tube, then add a small amount of 0.1M H<sub>2</sub>SO<sub>4</sub>. Add 6% H<sub>2</sub>O<sub>2</sub> solution to the test tube and observe what happens.

## Metal Sulfide


3. Prepare 3 test tubes with 1.0mL of H<sub>2</sub>S solution (0.2M), then add 0.5mL of 0.2 M CuSO<sub>4</sub>, Hg(NO<sub>3</sub>)<sub>2</sub>, and SnCl<sub>2</sub> to respective tubes.
4. Centrifuge the tubes, collect the precipitates, and observe their colors.
5. Add a small amount of 6M HCl to the CuS to see if the precipitate dissolves. After centrifuging and removing the supernatant, add a small amount of HNO<sub>3</sub> to the tube and the precipitate dissolving.
6. Add a small amount of sulfur powder to a test tube, followed by a small amount of 0.5M Na<sub>2</sub>S.
7. Heat the sample in a 60°C water bath and observe the color change. Add this solution into the SnS<sub>2</sub> precipitate test tube and observe the dissolution of the precipitate. Add 2.0 M HCl into the transparent solution and observe what happens.

# Procedure (Continued)

## Properties of $\text{Na}_2\text{S}_2\text{O}_3$


7. Add several drops of 0.5 M  $\text{Na}_2\text{S}_2\text{O}_3$  solution into a test tube containing 0.5 mL of 0.1 M  $\text{AgNO}_3$ . Keep adding drops until 5.0 mL of  $\text{Na}_2\text{S}_2\text{O}_3$  has been added. Observe the reaction.
8. Separately, add 5.0 mL of 0.5 M  $\text{Na}_2\text{S}_2\text{O}_3$  to a test tube. Add 0.5 mL of 0.1 M  $\text{AgNO}_3$  solution. Observe the reaction and compare to the previous reaction.

# Flow Chart



# Part 1

# Oxidative and Reductive activity of $\text{H}_2\text{O}_2$



1. Mix equal volume of 0.2 M  $\text{Pb}(\text{NO}_3)_2$  and  $\text{H}_2\text{S}$  solution.
2. Remove supernatant and add 30%  $\text{H}_2\text{O}_2$  solution.

The diagram illustrates the experimental procedure in three stages:

- Stage 1:** On the left, a green gas cylinder labeled "TYHJ H<sub>2</sub>S" and a brown glass bottle labeled "Lead (II) nitrate >99% CAS No. 10099-74-8 100G" are shown. Below them is a box containing the text "Prepared 0.2 M solutions of each".
- Stage 2:** A large grey arrow labeled "Equal volume" points from the first stage to a test tube containing a clear, colorless liquid.
- Stage 3:** A second large grey arrow labeled "Remove supernatant, add 30% H<sub>2</sub>O<sub>2</sub>" points from the test tube to a beaker. The beaker contains a dark, black precipitate at the bottom, with a clear liquid supernatant above it.

# Oxidative and reductive activity of $\text{H}_2\text{O}_2$ (cont.)

1. Add several drops of 0.01M  $\text{KMnO}_4$  in a test tube, followed by small amount of 0.1 M  $\text{H}_2\text{S}$ .
2. Add 6%  $\text{H}_2\text{O}_2$  solution in the test tube and observe.

Several drops 0.01M

Small amount of 0.1 M  $\text{H}_2\text{S}$




Add 6%  $\text{H}_2\text{O}_2$




# Part 2






# Metal sulfides



1. Respectively add 0.5 mL 0.2M  $\text{CuSO}_4$ ,  $\text{Hg}(\text{NO}_3)_2$ , and  $\text{SnCl}_2$  solution into three test tubes containing 1.0 mL  $\text{H}_2\text{S}$  solution (0.2 M).
2. Centrifuge, collect precipitations, and observe color difference.




## Metal sulfides (cont.)



3. Add small amount of 6M HCl to CuS test tube and see if precipitation will dissolve.
4. Add small amount of  $\text{HNO}_3$  into the test tube and observe precipitation dissolution.



Copper (II)  
sulfide





Small amount  
6M HCl



Small amount  
 $\text{HNO}_3$



## Metal sulfides (cont.)




1. Add small amount of sulfur powder into a test tube, followed by small amount of 0.5 M  $\text{Na}_2\text{S}$ . Heat sample in  $60^\circ\text{C}$  water bath and observe color change.




Sulfur powder


Small amount  
sodium sulfide



$60^\circ\text{C}$  water bath



## Metal Sulfides (cont.)



2. Add this solution into  $\text{SnS}_2$  precipitation test tube and observe precipitation dissolution. Add 2.0 M HCl into above transparent solution and observe results.

Tin (II) sulfide  
(in test tube)

2.0 M HCl

# Part 3

# Properties of $\text{Na}_2\text{S}_2\text{O}_3$

The diagram illustrates the reaction of sodium tetrathionate ( $\text{Na}_2\text{S}_2\text{O}_3$ ) with silver nitrate ( $\text{AgNO}_3$ ) in four stages, each shown in a test tube:

- Test Tube 1:** Contains 0.5 mL of 0.1 M  $\text{AgNO}_3$ .
- Test Tube 2:** After adding several drops of 0.5 M  $\text{Na}_2\text{S}_2\text{O}_3$ , a white precipitate ( $\text{Ag}_2\text{S}_2\text{O}_3$ ) forms.
- Test Tube 3:** After slowly adding a full 5 mL of 0.5 M  $\text{Na}_2\text{S}_2\text{O}_3$ , the white precipitate is more pronounced.
- Test Tube 4:** The color of the solution has changed to brown.

**Chemical Reaction:**



$$\text{Na}_2\text{S}_2\text{O}_3 + 2\text{AgNO}_3 \rightarrow \text{Ag}_2\text{S}_2\text{O}_3 + 2\text{NaNO}_3$$
$$\text{Ag}_2\text{S}_2\text{O}_3 + \text{H}_2\text{O} \rightarrow \text{AgS}_{(s)} + \text{H}_2\text{SO}_4$$

**Observations:**

- White ppt forms ( $\text{Ag}_2\text{S}_2\text{O}_3$ )
- Color changed to brown

# Properties of $\text{Na}_2\text{S}_2\text{O}_3$

New test tube



# Conclusion

**Accomplished:** While completing this lab we were able to illustrate the oxidative and reductive properties of  $\text{H}_2\text{O}_2$  and  $\text{H}_2\text{S}$ . We demonstrated the formation of multi sulfurs compounds. The reagents used in this experiment were  $\text{H}_2\text{O}_2$ ,  $\text{Pb}(\text{NO}_3)_2$ ,  $\text{H}_2\text{S}$ ,  $\text{KMnO}_4$ ,  $\text{H}_2\text{SO}_4$ ,  $\text{CuSO}_4$ ,  $\text{Hg}(\text{NO}_3)_2$ ,  $\text{SnCl}_2$ ,  $\text{HCl}$ ,  $\text{CuS}$ ,  $\text{HNO}_3$ ,  $\text{Na}_2\text{S}$ ,  $\text{SnS}_2$ ,  $\text{Na}_2\text{S}_2\text{O}_3$ ,  $\text{AgNO}_3$ . The techniques used were centrifuging, mixing, and observation.

**What was learned:** We learned that when  $\text{Pb}(\text{NO}_3)_2$  and  $\text{H}_2\text{S}$  are added together they form a black precipitate. When  $\text{H}_2\text{O}_2$  is added to this, the precipitate turns white. When  $\text{KMnO}_4$  and  $\text{H}_2\text{SO}_4$  are mixed, the solution is purple. When  $\text{H}_2\text{O}_2$  is added to this, the solution turns clear.  $\text{CuSO}_4$  starts as blue and when it is added to  $\text{H}_2\text{S}$ , it turns brown with a dark precipitate.  $\text{SnCl}_2$  is white/milky and when it is added to  $\text{H}_2\text{S}$ , it turns brown with precipitate.  $\text{CuS}$  and  $\text{HCl}$  together don't have a reaction. When  $\text{HNO}_3$  is added to this, the solution turns yellow. Sulfur powder with  $\text{Na}_2\text{S}$  in heat caused a color change from clear solution with white powder on top to a yellow solution with yellow precipitate on top. When this is added to  $\text{SnS}_2$ , the precipitate dissolves and the solution turns green. The addition of  $\text{HCl}$  makes the solution turn yellow with an orange precipitate.  $\text{Na}_2\text{S}_2\text{O}_3$  added to  $\text{AgNO}_3$  turns from clear to pale yellow with a white precipitate with the addition of addition  $\text{Na}_2\text{S}_2\text{O}_3$ . With more, the precipitate turns black and then dissolves.  $\text{AgNO}_3$  added to  $\text{NaS}_2\text{O}_3$  changes from a yellow precipitate to a brown precipitate.

**Any issues during the experiment:** Since the experiment was done virtually without calculations there were no issues, however a typical issue that may have occurred in the lab could have been heating and mixing for correct amount of time. If either of these were to be executed poorly then the reaction may not fully develop the precipitates and color changes meant to take place.

**Future recommendations/practical application:** Understanding the importance of oxidation/reduction reactions is essential especially when it comes to applying this knowledge to other research. Hydrogen peroxide is one of the most powerful oxidizers, but it can be used as a reducing agent as well.



# Postlab Questions

**1. What is the most advantage of using  $\text{H}_2\text{O}_2$  as an oxidant?**

$\text{H}_2\text{O}_2$  is advantageous to use as an oxidizing reagent due to its “purity” because when it breaks down it forms oxygen and water.

**2. If  $\text{H}_2\text{S}$  or  $\text{Na}_2\text{S}$  solution is placed for a long time, what may happen to the solution?**

The concentration will decrease as  $\text{H}_2\text{S}$  or  $\text{Na}_2\text{S}$  molecules bubble out of the solution.