

# Experiment 4: Preparation of Chloropentaamminecobalt (III) chloride and Nitropentaamminecobalt(III) chloride

Jenna Pralat, Emily McGuire, Kaylyn Spotton, Jessica Persaud, and Max Shi

We pledge our honor we have abided by the Stevens Honor Code.

# Procedure

## Preparation of Chloropentaamminecobalt (III) chloride

1. Combine 5.0 g of ammonium chloride and 30 mL of conc. aqueous ammonia in a 250 mL Erlenmeyer flask
2. Add stir bar to flask and start stirring
3. Slowly add 10 g of fine powder cobalt (II) chloride 6-hydrate.
4. While stirring, dropwise add 8mL of 30%  $\text{H}_2\text{O}_2$
5. Once reaction ceases, slowly add 30 mL of conc. HCl
6. Turn on hot plate to 85C, heat for 20 minutes
7. Cool to room temperature and Filter  $[\text{Co}(\text{NH}_3)_5\text{Cl}]\text{Cl}_2$
8. Wash crystals with ice water several times(max 20 mL), wash with 20 mL cold HCl (6 M)
9. Dry in 100 C oven for 2 hrs.

## Preparation of Pentaamminenitritocobalt(III) chloride

1. Prepare solution of 8 mL ammonia (aq) in 80 mL of water
2. Heat on hot plate, and a stir bar
3. Add 5.0 g of  $[\text{Co}(\text{NH}_3)_5\text{Cl}]\text{Cl}_2$  (dry)
4. Continue heating / stirring until colored product dissolves
5. Filter off dark brown/black cobalt oxide

# Procedure cont.


## Preparation of Pentaamminenitritocobalt (III) chloride

6. Cool solution to  $\sim 10^\circ\text{C}$
7. Add 2M HCl slowly, until neutral
8. Add 5.0 g of  $\text{NaNO}_2$  and 5 mL of 6M HCl.
9. Place in ice bath for 1 hr.
10. Filter precipitated  $[\text{Co}(\text{NH}_3)_5\text{ONO}]\text{Cl}_2$
11. Wash product with 25 mL ice water, 25 mL of alcohol, then dry on lab bench for 1 hr

## Preparation of Pentaamminenitrocobalt (III) chloride

1. Boil 20 mL of water
2. Add a few drops of ammonia (aq) and 2.0 g of  $[\text{Co}(\text{NH}_3)_5\text{ONO}]\text{Cl}_2$
3. As solution cools, add 20 mL of concd/ HCl
4. Filter crystalized  $[\text{Co}(\text{NH}_3)_5\text{NO}_2]\text{Cl}_2$  with Buchner funnel
5. Wash crystals with 13 mL alcohol
6. Allow to dry in air for 2 hrs


# Flowchart



ppt - precipitate

# 1. Preparation of pentaamminechlorocobalt (III) chloride

1. Make a solution of 5.0 g ammonium chloride in 30 mL concentrated aqueous ammonia in 250 mL Erlenmeyer flask.
2. Place flask on magnetic stirrer hot plate, turn on stirring, and add 10 g powdered cobalt (II) chloride hexahydrate.
3. With continued stirring, add 8 mL 30% hydrogen peroxide dropwise.



# 1. Preparation of pentaamminechlorocobalt (III) chloride

4. When evidence of reaction has ceased, slowly add 30 mL concentrated HCl.
5. With continued stirring, turn hot plate and adjust temperature to 85°C. Heat for 20 minutes at this temperature.
6. Cool mixture to room temperature.




Add 30 mL  
conc. HCl




85°C, 20 min,  
with stirring

Cool to  
room temp



Maroon color




Vacuum filter  
precipitate



When no more color  
change or gas evolution

# 1. Preparation of pentaamminechlorocobalt (III) chloride


7. Wash purple crystals with ice water not exceeding 20 mL, then 20 mL cold 6M HCl, and dry in 100°C oven for two hours.





## 2. Preparation of Pentaamminenitritocobalt (III) chloride



80 mL  
water



8 mL  
ammonium  
(aq)



heat a  
solution on  
the stirrer-hot  
plate





## 2. Preparation of Pentaamminenitritocobalt (III) chloride




Add 6 g of  
 $[\text{Co}(\text{NH}_3)_5\text{Cl}]\text{Cl}_2$   
while heating  
and stirring the  
solution



Heat until completely  
dissolved


## 2. Preparation of Pentaamminenitritocobalt (III) chloride

If a dark brown to black precipitate of cobalt oxide forms, filter it off



## 2. Preparation of Pentaamminenitritocobalt (III) chloride

Cool the filtrate which should be a clear solution to about  $10^{\circ}\text{C}$ .



Add 2 M HCl slowly while keeping the solution cold until it is just neutral to litmus.

## 2. Preparation of Pentaamminenitritocobalt (III) chloride

Should see salmon pink  
 $[\text{Co}(\text{NH}_3)_5\text{ONO}]\text{Cl}_2$


5.0 g of sodium nitrite

5 mL of 6M HCl

Add 5.0 g of sodium nitrite followed by 5 mL of 6 M HCl.

The diagram illustrates the preparation of Pentaamminenitritocobalt (III) chloride. It shows a beaker containing ice, with 5.0 g of sodium nitrite and 5 mL of 6M HCl being added. The resulting solution is salmon pink, indicating the formation of the product,  $[\text{Co}(\text{NH}_3)_5\text{ONO}]\text{Cl}_2$ .

## 2. Preparation of Pentaamminenitritocobalt (III) chloride



### 3. Preparation of Pentaamminenitrocobalt (III) chloride

Add few drops of aqueous ammonia

20 mL of water boiled


Measure and add 2.0 g of  $[\text{Co}(\text{NH}_3)_5\text{ONO}]\text{Cl}_2$

Add 20 mL HCl

Let solution cool

$[\text{Co}(\text{NH}_3)_5\text{ONO}]^{2+} \xrightarrow{\text{heat}} [\text{Co}(\text{NH}_3)_5\text{NO}_2]^{2+}$

### 3. Preparation of Pentaamminenitrocobalt (III) chloride



# Conclusion

Accomplished: In this experiment, pentaamminechlorocobalt (III) chloride was prepared using  $\text{H}_2\text{O}_2$  to oxidize  $\text{Co}^{2+}$  into  $\text{Co}^{3+}$  and  $\text{HCl}$  for the ligand exchange. The pentaamminechlorocobalt (III) chloride was then used to prepare pentaamminenitritocobalt (III) chloride using  $\text{NaNO}_2$  for a ligand exchange. Lastly, the pentaamminenitritocobalt (III) chloride was isomerized into pentaamminenitrocobalt (III) chloride by heating, since these two compounds are linkage isomers.

Learned: Throughout this lab, we learned about complex ions and coordination compounds. The different types of structural isomers (coordination and linkage) were learned through the experiment, as well as how they can be prepared in the laboratory through ligand exchange reactions and isomerization.

Any Issues: The main issues in the experiment were not obtaining enough precipitate of the coordination compounds to proceed to the next step in the procedure due to inefficient reactions or improper heating during the reactions.

Future Recommendations / Applications: To increase the amount of product obtained in each step of the reaction the temperature should be very closely monitored throughout the process and adjusted to be at the optimal level for each step.



# Post Lab Questions

1. What is the name of the following type of reaction?
  - a.  $[\text{Co}(\text{NH}_3)_5\text{H}_2\text{O}]^{3+} + 3\text{Cl}^- \rightarrow [\text{Co}(\text{NH}_3)_5\text{Cl}]\text{Cl}_2 + \text{H}_2\text{O}$
  - b. This reaction is a ligand exchange reaction where the counterion and ligand switch, forming the coordination isomer.
2. Why did ammonia not get displaced in the above reaction?
  - a. Ammonia is not displaced because its interaction with cobalt is more stable than that of water, so the water can be easily replaced, whereas the ammonia can not.
3. Name two different techniques which can be used to identify whether it is nitrito complex  $[\text{Co}(\text{NH}_3)_5\text{ONO}]^{2+}$  or nitro complex  $[\text{Co}(\text{NH}_3)_5\text{NO}_2]^{2+}$ ?
  - a. The first technique is by observation of the different colors. The nitrito complex is a red color and the nitro complex is a yellow color. The second technique is by melting point determination, the H-bonds that are possible in nitro will increase the melting point because the solid nitro complex is more stable than the nitrito complex.